




Premium Control hub. Your premium needle.



TSK Laboratory Europe BV
De Scheper 309
5688HP Oirschot
The Netherlands

europa@tsklab.com
www.tsklab.com
+31 85 760 7630

Premium injection
solutions for
aesthetics

Up to 50% less drag force.

The PRC Control Hub delivers refined precision and enhanced control for dermal filler injections. Its ergonomic and stable hub design improves handling, enabling practitioners to administer treatments with confidence and consistency.

At its core is TSKiD™ Large technology, which features an enlarged inner diameter that reduces extrusion force and supports smoother, more controlled injections into high-


viscosity gels. This optimized needle flow significantly improves injection comfort and accuracy.

Mechanical security is reinforced by the hub's construction from an extremely hard polymer, which prevents flexing under high pressure.

This rigidity minimizes the risk of leakage and pop-off, ensuring dependable performance throughout demanding procedures.

The PRC Control Hub's traditional hub geometry maintains full compatibility with standard syringes, supporting seamless integration into existing practice workflows without compromise.

Designed for precision and reliability, the PRC Control Hub enhances stability and control in advanced aesthetic applications, helping clinicians deliver consistent, high-quality results.






Lancet bevel cut

The ultra sharp lancet bevel cut makes it easier to insert the needle.



Thin wall barrel

Enlarged inner diameters ensure an enhanced flow rate and less extrusion force.



Hard polymer hub

The hub is made of polycarbonate. A harder polymer making the hub stronger.



Reliable sharpness.



Optimal patient comfort.

The innovative TSK STERiJECT™ bevel cut is known for its exceptional sharpness, providing a precise, efficient and comfortable injection with each use.

Sterile packaging.



Safe hard case capsule.

TSK's safety capsule for hypodermic needles, allows for easy and safe needle-to-syringe assembly with minimal risk of needle stick injury.

Available sizes for the Premium Control needle:

Product code	TSKiD Inner Diameter	Gauge	Length
PRC-30013I-100	Large	30G	13mm (1/2 inch)
PRC-27038I-100	Large	27G	38mm (1 1/2 inch)
PRC-27019I-100	Large	27G	19mm (3/4 inch)
PRC-27013I-100	Large	27G	13mm (1/2 inch)
PRC-25013I-100	Large	25G	13mm (1/2 inch)

Premium injection solutions for aesthetics

TSK Laboratory Europe BV
De Scheper 309
5688HP Oirschot
The Netherlands

Order now
www.tsklab.com

